

MICROSCOPE COMPARISON CHART

BEST SELLERS

SPECIALIZED



	GOOD	BETTER	BEST				
MODEL	Revelation III	i4 Infinity	InnoVation	InnoVation with BioVIEW	i4 Semen Evaluation	Mi5 Dual Head	Z4 Embryo-GLO Zoom Stereo
MODEL NUMBER	J0334AB	J0334Q	J0334iB4	J0334iT4X	J0334QSE	J0334MI5DH	J0755Z4BEG
HEAD	Binoc (Trinoc available)	Binoc (Trinoc available)	Binoc (Trinoc available)	Trinoc	Binoc (Trinoc available)	Trinoc + Binoc	Binoc (Trinoc available)
STAGE	Mechanical	Mechanical	Rackless Mechanical	Rackless Mechanical	Mechanical - Heated	Round, Rotates 360, Mechanical Optional	Frosted Stage Plate
OBJECTIVES	DIN Achromatic 4-10-40-100 oil	Infinity S-PLAN 4-10-40-100 oil	Infinity PLAN 4-10-40-100 oil	Infinity PLAN 4-10-40-100 oil	Infinity PLAN PHASE 10-20-40-100 oil	Infinity PLAN 4-10-20-40-100 oil	Zoom 7x to 45x
NOSEPIECE	Quadruple Forward-Facing	Reverse Quadruple	Reverse Quadruple	Reverse Quadruple	Reverse Quadruple	Reverse Quintuple	N/A
FEATURES	Affordable, Durable, LED Light	High Res, 95% Flat-field, LED Light	High res, flat-field, lab-grade, rackless stage, LED light	HD Camera and HD Monitor mounted onto Trinoc InnoVation	Heated stage, Phase for semen, also brightfield	Two people viewing same sample	Zoom optics, tilt-mirror for contrast, higher mag options
APPLICATIONS	Urine, Fecal, Blood, Basic Cytology	Urine, Fecal, Blood, Routine Cytology	Urine, Fecal, Blood, Heavy Cytology	Client Education, Tech Training, University, Email To Pathologist	Semen Evaluation, Urine, Fecal, Blood, Cytology	University, In-Clinic Tech Training, Dermatology MOHS Technique	Embryo Transfer, Dissections

*5 Objective Model Available

*Binoc + Binoc Model Available