

OLYMPIC VERSA-LIFT™

One electric lift that does 3 jobs



GURNEY

Lightweight, easy to roll and steer

TREATMENT

IV Pole and Instrument Tray fit in corner sockets



CURBSIDE

Remove patients from cars easier, safer with electric lift and lift-off stretchers.



Hands-on care

Raise lift for transport with assistant.

Hidden handle
Slides out. One person can safely move patients outdoors.



New Model Versa-Lift
with
QUICK-LOCK
Casters



One lever locks/unlocks all 4 casters

Saves time No more locking and unlocking one caster at a time.

Safer Locks all 4 casters securely. No chance of an unlocked caster. 300 lb. patient capacity.

2 Brake Levers One on each side of the Versa-Lift.

2 Foot Switches For electric lift, one on each side of Versa-Lift.

Instrument Tray Included on swing arm, drops into sockets.

Order Model 50705

More >>>

OLYMPIC VERSA-LIFT™

For lifting and transporting Saves time Improves patient care



Features

Protects against back injury Ergonomic electric lift. No more lifting or carrying heavy patients. Helps your hospital meet safety standards. Maximum height 40" for standing comfort.

Saves time Takes fewer people and less time to transport heavy patients.

Improves patient care A gurney that combines an electric lift with lift-off stretchers. Less traumatic for patients than awkward hand carrying, especially for injured or compromised patients.

Safety compliant Helps your hospital meet current safety standards.

Rolls easily Lightweight, only 68 pounds, half the weight of some gurneys. Large, locking casters. (See reverse side for new Versa-Lift with QUICK-LOCK Casters.)

Turns on a dime Steering handles at both ends. Maneuvers in tight spaces.

Low-boy handle With lift in lowered position, a hidden handle slides out so one person can safely transport a patient in the low-boy position if necessary.

Outdoor use Makes removing patients from cars easier, safer with electric lift and lift-off stretchers.

Quick recharge battery With charging light.

Versatile 3 products for price of one – a gurney, treatment table, and curbside lift.

Watch video demonstration www.OlyVet.com 

Lift-off Stretchers Ergonomic and safe



The top of the Versa-Lift is a stretcher which is held securely on 4 posts but easily lifts off for transferring patients to x-ray and other tables. Saves back strain and complies with safety recommendations. Easier, quicker transfers. Safer for patients and staff. Two stretchers available.



Hard Stretcher

For treatment and exams. Molded of rigid reinforced polymer. Durable. Wipes clean. Non-slip surface. 2 Velcro straps.



Soft Stretcher

Reinforced fabric is medical-grade, anti-microbial Herculite®. Wipes clean. X-ray transparent. Strong. 2 Velcro straps.

Specifications

Dimensions 52" x 24".

Height range 13" to 40".

Casters Large 5" casters with new Quick-Lock central locking system, or individual locks.

Net weight 68 pounds.

Steering Handles At both ends.

IV Pole Included all models.

Instrument Tray on swing arm included on Model 50705.

Foot Switches 2 foot switches on Model 50705. One switch on other models.

Electrical 120v/60Hz.

Order from your Distributor



Olympic Versa-Lift with Quick-Lock Casters. Complete with Hard Stretcher, Soft Stretcher, IV Pole, Instrument Tray, and 2 Foot Switches **Model 50705**

Olympic Versa-Lift with Hard Stretcher, Soft Stretcher, and IV Pole **Model 50701**

Olympic Versa-Lift with IV Pole and Hard Stretcher only.... **Model 50703**



Instrument Tray Accessory fits on Versa-Lift **Model 50780**

3-Year Olympic Warranty Includes all parts and labor