



SIMPLIFYING YOUR CONTROLLED SUBSTANCE INVENTORY MANAGEMENT PROCESS!

The CSI-360 combines a two-part authentication process, novel software and RFID tracking technology to provide secure access and 100% accurate recording of all controlled substance use in the veterinary practice. All records are captured electronically for a precise account of all cabinet access and activity, saving time and eliminating data entry errors. Detailed inventory reports are easily retrievable and readily available to ensure DEA compliance.

Improve the Profitability of Your Veterinary Practice

- Save time and costs associated with manual logging and recordkeeping
- Use your existing DEA-certified cabinet, safe, fridge or freezer
- Save costs by optimizing the appropriate inventory needed on hand

Enhance Security to Protect Your Practice

- Secure access with electronic lock and dual authentication to prevent unauthorized access
- Complete accountability of all cabinet activity, including shortages and misuse
- Accurate recordkeeping and audit trails to prevent theft, ensure compliance and avoid fines

Increase Efficiency and Productivity

- Reduce manual labor with automatic and precise tracking and recording
- Eliminate data entry errors with paperless recording
- Fast and easy access to inventory reports



The UID narcotics cabinet system was an affordable game changer for our practice. The CSI-360 modernized our operations by increasing security and reducing the overall liability for our doctors."



The CSI-360 provides enhanced security features and paperless recording for efficient inventory management of controlled substances and other medications in the veterinary practice.



Unique Capabilities

- ✔ Secure access with two-factor verification
- ✔ Authorize or restrict user access
- ✔ Accurate inventory of all contents
- ✔ Paperless recording for accurate documentation
- ✔ Automatic and electronic tracking from acquisition to administration
- ✔ Detailed reports: inventory, activity, user access and more
- ✔ Network Multiple Cabinets
- ✔ All documents are readily available for DEA audits
- ✔ Compatible with all cabinets, safes, refrigerators, and freezers from most manufacturers



System Specifications

- Cabinet/Safe (optional): use your existing DEA-certified cabinet/safe
- Electronic lock: software-controlled, secondary lock for extra security
- Computer: Windows surface PC with SQL database inventory management software
- Key card, fob or wristband for fast and touch-free user access
- Scanner: RFID reader for user access and item identification
- Digital scale: for accurate tracking by weight and density
- Labels: RFID labels for unique identification of stored items



Flexible Purchasing Options to Accommodate Your Practice's Budget

CONTACT US FOR DETAILS 

