


Kanine and Kitty Kollar®

The collar provides a simple, practical solution for everyday management of esophagostomy (e-tube) small animal patients. It fits securely around the neck, features a dedicated tube opening, and includes a Velcro™ strap to secure the tube when not in use. Supplied circular cotton pads protect the wound, absorb drainage, and shield against contamination. This system eliminates the need for owner-applied bandaging and reduces repeat visits for re-bandaging.

- Machine washable collars provide superior cleanliness of site
- Ease of application and removal equals convenience
- Lowered risk of tube removal by animal
- Decreased risk of wound trauma (animals scratch at site less; less chance that owners will bandage too tightly)
- Can leave in place for radiographs, MRI, chemotherapy, other procedures
- Send home a more professional and compassionate presentation and aesthetically pleasing
- Each collar is packaged in a plastic bag with a laminated information card and 6 protector pads (a one-to two week supply)



10% off
March 1 -
April 30,
2026!

Item No.	Size	Breed	Weight	Neck Size	Length (collar + tab)	Width	End of Collar to Tube Hole	Top Edge of Collar to Tube Hole
350.79012.2	X-Small	Kittens/Toy Breeds (Chihuahua, Yorkie)	4-7 lbs	4"-7"	5.5" + 2"	2"	1.5"	3/4"
350.79006.2	Petite	Sm Cats/Dogs (Yorkie, Poodle)	6-9 lbs	8"-9"	9" + 2"	2"	1.5"	3/4"
350.79007.2	Small	Cat/Small Dog (Lhasa, Pekingese, small terriers)	10-16 lbs	10"-12"	11" + 2"	2"	1.5"	3/4"

Notes:

- Long-haired breeds may need to go up a size
- More sizes available for larger dogs
- Color/pattern is variable

