



**YOUR PARTNER IN
VETERINARY PRODUCTS.**

EXAM GLOVES



Latex: Exam grade, PF latex, high level of touch sensitivity, wear for extended periods of time, very elastic and strong, biodegradable, latex gloves are a natural material, made out of rubber. High level of barrier protection.

Vinyl: PF Vinyl, latex-free, looser fit, short term option, lightweight alternative, great for everyday use in many aspects, general care to be disposed of after each use.

Nitrile: Exam grade, PF nitrile, latex-free, textured finger tips, strong, beaded cuff, snug and comfortable fit, high level of sensitivity, wear for extended periods of time, made out of synthetic rubber, and are a good alternative when latex allergies are of concern. High level of barrier protection.

All exam gloves are powder free.

Size	Palm Width (in)	Hand Height
S	2.75" +/- 0.4"	9.5"
M	3.10" +/- 0.4"	9.5"
L	3.75" +/- 0.4"	9.5"
XL	4.30" +/- 0.4"	9.5"



LATEX:

ITEM	SIZE	COUNT
000.30040.2	XL	100 GLOVES
000.30042.2	LG	100 GLOVES
000.30044.2	MED	100 GLOVES
000.30046.2	SM	100 GLOVES

VINYL:

ITEM	SIZE	COUNT
000.30050.2	XL	100 GLOVES
000.30052.2	LG	100 GLOVES
000.30054.2	MED	100 GLOVES
000.30056.2	SM	100 GLOVES

NITRILE:

ITEM	SIZE	COUNT
000.30070.2	XL	200 GLOVES
000.30072.2	LG	200 GLOVES
000.30074.2	MED	200 GLOVES
000.30076.2	SM	200 GLOVES

For our full product line, or to place an order,
visit midwestvetsupply.com or call 800.MIDWEST (643-9378)

Distributed Exclusively by:

